

Sell Sheet | On-chip Sort

Value Proposition: The On-chip Sort provides damage / stress free sorting on a disposable, easy to use microfluidic chip, with high purity and recovery. This provides a unique benefit to those researchers with rare or fragile cells that need precision sorting that won't negatively impact their cell samples.

Flow cytometry is a technique used to detect and measure physical and chemical characteristics of a population of cells or particles

Major Features of On-chip Sort:

Damage free sorting

- Low flow pressure (0.3 – 3 psi)
- Gentle air pulses for cell deflection
- No electrostatic charge
- Free to use any sheath fluid (culture media, any kind of buffer, oil, etc.)
 - On-chip does offer a proprietary sheath fluid that is well tested and characterized for those customers that don't require very specialized sheath fluid

High recovery (>90%)

- This is how much of the target cell is being collected from the starting sample.
- Ex: If our starting sample has 100,000 cancer cells, we're interested in how much of that 100,000 is gathered after sorting. In this case, it should ideally be >90,000 cells

High purity (>95%) depending on cell concentration

- This is the percentage of sorted cells that were the target cell of interest and not miscellaneous cells. In other words, how accurate the sort was.
- In the previous example, if we managed to sort 90,000 cells, we should only have an extra 2,000-3,000 non-target cells that were sorted alongside our target of interest.

Lasers:

- Up to three types including **Blue** (488 nm) as standard and two selected from **Violet** (405 nm), **Green** (561 nm), or **Red** (638 nm) can be used (additional types are also available upon request).

→ Target Market and their needs

- **Flow cytometry core labs**
- **3D cell culture researchers**
 - 3D cell aggregates (e.g. spheroids/organoids) are traditionally difficult to sort with conventional cell sorters due to their size.
 - It is possible to sort spheroids/organoids in the 150 µm channel microfluidic chips
- **Microbiologists**
- **Researchers with very rare or very fragile cells**
 - Samples with target cells that are <100:100,000
 - Cells that might differentiate or change gene expression under stress
 - Rare clinical samples
- **Users performing droplet sorting**
 - water-in-oil droplets
 - Gel droplets



On-chip Sort

Sell Sheet | On-chip Sort

Value Proposition: The On-chip Sort provides damage free/stress free sorting on a disposable, easy to use microfluidic chip, with high purity and recovery. This provides a unique benefit to those researchers with rare or fragile cells that need precision sorting that won't negatively impact their cell samples.

Additional Major Features of On-chip Sort:

Microfluidic chip


- Disposable
- All sorting activity occurs on the chip
 - No clean up
 - No cross contamination
 - Easy to run subsequent sorting runs on the target sample or waste sample
- High clarity plastic for microscopy
- 80µm and 150µm channel sizes for a variety of samples – including large organoids and gel droplets

Ease of use

- Very little maintenance or labor when using the hardware or inserting the microfluidic chip
- Straightforward software
- Because of the easy-to-use hardware and software, multiple users can be quickly taught to use the On-chip sort
- **Start-up takes as little as 5 minutes; no calibration is required.**

Application - cell sorting and flowcytometry

- Rare or fragile cell isolation (**key feature**)
- Large cells, organoids/spheroids, tissue samples (**key feature**)
 - Up to 50µm targets in the 80um chip
 - Up to 120µm targets in the 150um chip
- Bacterial sorting
- Gel and oil encapsulation – cell(s) **encapsulated** in droplets, usually remaining viable and active within the droplet microenvironment (**key feature**)
 - Gel droplet, Nano vials
 - Water-in-oil droplets



Microfluidic chip

→ Target Competitors

- **BD Biosciences – FACSaria:**
These units have higher throughput. A common conventional cell sorter
- **Miltenyi – MACSQuant:**
These units have higher throughput
- **Sony – FX500, SH800S, MA900:**
Some units use a different type of microfluidic chip-based sorting technology. Some units are heavily focused on GMP
- **Namocell - Pala:**
Similar microfluidic chip technology. Less features, but cheaper list price
- **Nanocellect - WOLF:**
Similar microfluidic chip technology. Less features, but cheaper list price