



Equip Your Laboratory With The
SAFEST & MOST PERSONALIZED GROSSING BENCH



www.mopec.com



(800) 362-8491



info@mopec.com

Mopec Maestro™ Grossing Station

Mopec understands that each pathology professional has their own workspace preferences, and each laboratory holds their own safety, ergonomic and workflow requirements. With the Maestro's distinctive Ensemble Pegboard™, dozens of accessories and flexible safety options, every station can be quickly customized to suit each individual.



Premium Worksurface



CUTTING & RINSING SURFACE

Each station features two removable panels with one perforated plate and one specimen cutting board with removable legs. The panels are capable of sliding across set rails for workspace flexibility and easy cleaning. The sink can also be adjusted for left or right handed users.

HEIGHT ADJUSTABILITY

Maestro workstations are ADA compliant, and they adjust to the tallest height in the industry.



Intuitive Controls & Device Management

EASY-TO-USE TOUCHSCREEN

With a quick tap, all Maestro features can be controlled using the new, integrated touchscreen. All user settings can be saved with individual user profiles to allow professionals to promptly adjust to their height, lights and ventilation preferences.

WIRE CONTAINMENT

Four outlets are provided on each bench with three USB compatible and one GFCI. Cable covers are incorporated into the peg board and on the side of each station, therefore eliminating the look of messy cords.



Top-Notch Ventilation & Safety



FRONT AIR SYSTEM

This patented on-demand technology pushes air from the front of the cutting surface through the back of the area, creating a laminar flow pattern. This constant ventilation protects the user from dangerous formalin exposure.

SAFETY NOTIFICATIONS

A series of audible beeps and on-screen warnings alert users of possible harmful conditions involving air flow, formalin collection and other hazards.



Harmonize Your Workspace with Optional Accessories

Pegboard Accessories

Movable options for precise workstation personalization

Add-on pegboard wings allow for more storage space.

- Organizer Bins
- Paper/File Holder
- Bar Code Holster
- Shelves
- Cell Phone Holder
- Magnetic Tool Holder
- Staining Dye Holder
- Specimen Scale Shelf
- Towel & Glove Box Holder
- SaniPath Wipe Dispenser
- Glasses Holder
- Face Shield Holster
- Cable Organizer
- Center Monitor Mount



Optional Drawer System

Sliding options for quick reach accessibility

- Side Drawer
- Middle Drawer
- Small Over-the-Sink Drawer



LED lightstrip for an illuminated worksurface.



Formalin Fume Detection

The Mopec Maestro™ is the only grossing bench to offer optional formalin detection. This advanced technology alerts the user when an adjustable exposure threshold is met or exceeded. This monitoring system is ideal for laboratories working with large specimens and sizable amounts of formalin to ensure ultimate user safety.

Other Accessories

Additional selections for an enhanced workspace

- Cassette Labeling Shelf
- Side Monitor Arm
- Spatter-Free Formalin Dispenser
- Ventilated Trash
- CPU Bracket
- Keyboard Arm
- Integrated Camera
- Integrated Power Strip
- Mouse Tray
- Safety Splash Shield
- LED Light Strip
- Magnifier
- Task Light
- Keyboard Tray

Mopec Maestro™ Stations are ideal for:

- Hospital Pathology
- Private Pathology
- Animal Research
- Universities
- Histology
- Government Research

Why choose the Mopec Maestro™?



FUTURE FLEXIBILITY

Should laboratory space or staffing needs change, Maestro™ benches can be altered within your facility. Mopec technicians have the ability to change the sink location, add peg board wings or easily replace outdated parts, among other tasks.

QUALITY CONSTRUCTION

Each order is thoroughly tested and rigorously inspected against the highest industry standards. Every unit is issued with a standard one year warranty.

KNOWLEDGEABLE REPRESENTATIVES

Mopec understands that every facility is unique. Our customer service and sales representative teams are educated in common facility challenges and can assist in customizing your equipment.

Made in the United States

Worldwide Shipping

Over 25 Years of Experience



www.mopec.com
(800) 362-8491
info@mopec.com