



Ventilated Tissue Storage

No More Loading and Unloading



Top Vented
Air Duct
Connection

Vinyl
Coated
Steel Racks
w/Adjustable
Shelves

Adjustable
Floor Section
w/Leveling
Feet

Heavy-Duty
5" Locking
Swivel
Casters

Mopec The Better Product People

Call Us at: **1-800-362-8491** or Visit Us on the Web at www.mopec.com

VENTED TISSUE STORAGE AND TRANSPORT

A safe and highly functional method of tissue storage and transportation, the Ventilated Tissue Storage System is the practical answer to your storage needs. The system includes a portable shelving rack, container trays, transport cover and ventilated cabinet. With a portable rolling rack available, it eliminates the need to load and unload trays. For convenience, small containers may be placed in container trays with formalin neutralizing pads for transport and containment. The Mopec Ventilated Tissue Storage System is easy to adapt to most laboratory ventilation systems with an easily accessible 6" diameter duct stub located on the top of the cabinet. For convenience and functionality, the door to the cabinet is clear heavy-duty Lexan and may be hinged left or right. Four leveling feet and an adjustable stainless steel floor section allow convenient and easy installation.

CABINET FEATURES:

- Top ventilation duct connects to laboratory ventilation
- Heavy-duty 1/2" clear Lexan door (hinged left or right)
- Four leveling feet with an adjustable floor section
- All stainless steel for easy cleaning

Product No. LD500

PORTABLE SHELVING FEATURES:

- Fully adjustable shelves (1.5" intervals) allow for virtually any size specimen container
- Stable, vinyl coated steel racks
- Removable container trays with formalin neutralizing pads
- Large 5" heavy-duty swivel casters (2 locking)

Product No. LD501



Discreet form-fitting cover for LD500,
order Product No. JZ011

"Everyone talks about time savings, I talk about how much better my back feels not having to load and unload racks all day!"
— Morgue Technician

For a complete catalog detailing our products and services, please call or visit our web site



21750 Coolidge Highway • Oak Park, MI 48237
Toll Free: 800-362-8491 • Tel: 248-291-2040 • Fax: 248-291-2050

www.mopec.com