



PH-LRP-HC-SSP-49FA

The Plus Series Stainless Steel Pharmacy Refrigerators and Freezers from Lab Research Products are purpose built for laboratory and medical applications.

All models feature a digital microprocessor temperature controller, a forced air directional refrigeration system, and an automatic cycle-defrost design that prevents product damaging temperature spikes.

These products are ETL listed and feature environmentally friendly hydrocarbon natural refrigerant, R290.



FEATURES

- Digital display microprocessor temperature controller
- Adjustable operating temperature:
Refrigerator: 36°F to 46°F
Freezer: 4°F to -13°F
- Temperature Monitor Device - complies with the current CDC guidelines, 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory, °F/°C switchable, field installable, and visual & audible temperature alarms
- Audible and visual Hi/Lo temperature alarm
- Remote alarm contacts
- Forced draft air circulation
- Keyed door lock
- Magnetic door gaskets for positive seal
- LED interior lights are safety shielded and switch controlled
- Probe access port (3/4”)
- Stainless steel exterior
- Hydrocarbon natural refrigerant (R290)
- Pharmacy refrigerator/freezer toolkit and temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- 2 year parts and labor warranty, plus an additional 3 year compressor parts warranty
- Cycle defrost
- 115V/60 Hz
- ETL listed



PH-LRP-HC-SSP-23G

SPECIFICATIONS

Refrigerator Models	Ft ³	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
PH-LRP-HC-SSP-23	23	Cycle	1 Solid	3 Adj.	26 7/8	34 3/4	83 1/4	R290	1/5	3
PH-LRP-HC-SSP-23G	23	Cycle	1 Glass	4 Adj.	26 7/8	34 3/4	83 1/4	R290	1/5	3
PH-LRP-HC-SSP-49	49	Cycle	2 Solid	6 Adj.	54	34 3/4	83 1/4	R290	1/4	4.5
PH-LRP-HC-SSP-49G	49	Cycle	2 Glass	8 Adj.	54	34 3/4	83 1/4	R290	1/4	4.5
PH-LRP-HC-SSP-72	72	Cycle	3 Solid	9 Adj.	81	34 3/4	83 1/4	R290	1/2	7.5
PH-LRP-HC-SSP-72G	72	Cycle	3 Glass	12 Adj.	81	34 3/4	83 1/4	R290	1/2	7.5
Freezer Models										
PH-LRP-HC-SSP-23FA	23	Auto	1 Solid	3 Adj.	26 7/8	34 3/4	83 1/4	R290	3/4	8
PH-LRP-HC-SSP-49FA	49	Auto	2 Solid	6 Adj.	54	34 3/4	83 1/4	R290	3/4	9



Lab Research Products Plus Series Stainless Steel Pharmacy Refrigerators and Freezers take a large step forward in performance from past general purpose units; offering a digital temperature display and alarm module coupled with a microprocessor temperature controller.

