


The Basic Series Manual Defrost Freezers feature a digital microprocessor temperature controller allowing for the precise temperature management necessary for critical samples and supplies. The temperature display and alarm module offers product security with a monitoring and alarm system. Real-time sample simulated temperatures are constantly displayed and supported by audible/visual high and low temperature alarms. Remote alarm dry contacts and an access port are provided to simplify connections to existing facility monitoring systems.

Manual defrost units have been designed to operate without product damaging auto defrost cycles. **Periodic manual defrosting is required.**



DAI-HC-MFP-14

DAI-HC-MFP-20






FEATURING

- Microprocessor temperature controller
- Adjustable temperature range: -15°C to -25°C
- Manual defrost
- Remote alarm contacts
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- Recessed door handle
- Probe access port
- ETL Listed
- Two year parts and labor warranty, plus an additional three year compressor parts warranty
- Models require 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure



- Digital temperature display
- Audible and visual temperature alarms
- Resettable Min/Max
- °C/°F switchable

| Freezer Model | Ft ³ | Set Point | Defrost | Doors | Shelves | Width" | Depth" | Height" |
|---|-----------------|-----------|---------|---------|---------|---------|---------|---------|
|  DAI-HC-MFP-14 | 14 | -20°C | Manual | 1 Solid | 5 Adj. | 23 5/8" | 24 3/8" | 73" |
|  DAI-HC-MFP-17 | 17 | -20°C | Manual | 1 Solid | 5 Adj. | 30 3/4" | 28 5/8" | 67 3/4" |
|  DAI-HC-MFP-20 | 20 | -20°C | Manual | 1 Solid | 5 Adj. | 30 3/4" | 29" | 74 5/8" |

Specifications are subject to change