

Plus Series Manual Defrost Freezers feature a microprocessor temperature controller, allowing for precise temperature control, recovery and maintenance versus mechanical style thermostats. Manual defrost freezers also feature insulating inner doors and a manual defrost operating system to limit product-damaging temperature spikes.

In addition to the elimination of damaging auto defrost cycles, the inner doors help reduce cold air loss during main door openings and greatly improve temperature recovery. **Periodic manual defrosting is required.**

FEATURING

- Microprocessor Temperature Controller
- Adjustable temperature range:
-20°C Models: -15°C to -25°C
-30°C Models: -27°C to -33°C
- Manual defrost
- Remote alarm contacts
- Inner "fast-freeze compartments" with doors for temperature maintenance and recovery (30 ft³ only)
- Keyed door lock
- Magnetic door gasket for positive seal
- Caster platform (4") (30 ft³)
- Casters (2 ¾") (15 ft³)
- Probe access port
- CFC-free refrigerant, foam insulation and packaging (30 ft³)
- UL/C-UL listed (30 ft³)
- Two year parts and labor warranty, plus an additional three year compressor parts warranty



DAI-MFP-3030




DAI-MFP-15-C



- Digital temperature display and alarm module with battery back-up to support temperature monitoring and alarms
- Sample simulated digital temperature display; .1° Resolution, °C/°F switchable
- Full suite of alarms including - audible and visual high/low temperature alarms, sensor error alarm, power failure alarm
- Resettable 1-point min/max history
- Alarm validation (audible, visual, and remote alarm)
- Alarm mute with ring back



Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates. (15ft³ only)

Freezer Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
DAI-MFP-30	30	-20°C	Manual	1 Solid	5 Insulated	35 ¾	34 ¼	82 ½*
DAI-MFP-3030	30	-30°C	Manual	1 Solid	5 Insulated	35 ¾	34 ¼	82 ½*
Chest Freezer Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 DAI-MFP-15C	15	-20°C	Manual	1 Hinge	2 Baskets	56	29 ½	35 ⅞†

Specifications are subject to change

*Height includes 4" caster platform

†Height includes 2 ¾" casters