

The Basic Series Glass and Solid Door Laboratory Refrigerators are purpose built for scientific and medical applications. To improve performance over past general purpose units, the Basic Series Glass & Solid Door Laboratory Refrigerators feature a digital microprocessor temperature controller allowing for the precise temperature management necessary for critical samples and supplies. A forced air directional refrigeration system provides superior temperature uniformity and recovery after door openings. An automatic cycle-defrost design prevents product damaging temperature spikes.

FEATURING

- Microprocessor temperature controller
- Adjustable temperature range: 1°C to 10°C
- Cycle defrost
- Remote alarm contacts
- Forced draft air circulation
- Keyed door lock
- Magnetic door gasket for positive seal
- LED interior lights
- Swivel casters (3 ½") (select models only)
- Probe access port
- ETL Listed
- One year parts and labor warranty, plus an additional four year compressor parts warranty
- Select models require 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure (14, 17, and 20 ft³ models only)



DAI-HC-RFP-14

DAI-HC-LB-47





DAI-HC-SLB-49



- Digital temperature display
- Audible and visual temperature alarms
- Resettable Min/Max
- °C/°F switchable



Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates.

Glass Door Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 DAI-HC-LB-12	12	4°C	Cycle	1 Swing Glass	3 Adj.	25	29 ¾	65 ¾*
 DAI-HC-LB-16	16	4°C	Cycle	1 Swing Glass	3 Adj.	25	29 ¾	79*
 DAI-HC-LB-23	23	4°C	Cycle	1 Swing Glass	4 Adj.	26 7/8	34 ¾	81 ¾*
 DAI-HC-LB-26	26	4°C	Cycle	1 Swing Glass	4 Adj.	28 ¾	36 ½	81 ¾*
 DAI-HC-LB-33	33	4°C	Cycle	2 Sliding Glass	8 Adj.	39 5/8	33 ¼	81 ¾*
 DAI-HC-LB-36	36	4°C	Cycle	2 Swing Glass	8 Adj.	39 5/8	34 ¾	81 ¾*
 DAI-HC-LB-47	47	4°C	Cycle	2 Sliding Glass	8 Adj.	54 1/8	33 ¼	81 ¾*
 DAI-HC-LB-49	49	4°C	Cycle	2 Swing Glass	8 Adj.	54	34 ¾	81 ¾*
 DAI-HC-LB-69	69	4°C	Cycle	3 Sliding Glass	12 Adj.	81	33 ¼	81 ¾*
 DAI-HC-LB-72	72	4°C	Cycle	3 Swing Glass	12 Adj.	81	34	81 ¾*
Solid Door Model								
 DAI-HC-SLB-12	12	4°C	Cycle	1 Swing Solid	3 Adj.	25	29 ¾	65 ¾*
 DAI-HC-RFP-14	14	4°C	Cycle	1 Swing Solid	5 Adj.	23 5/8	24 3/8	73
 DAI-HC-SLB-16	16	4°C	Cycle	1 Swing Solid	3 Adj.	25	29 ¾	79*
 DAI-HC-RFP-17	17	4°C	Cycle	1 Swing Solid	5 Adj.	30 ¾	28 5/8	67 ¾
 DAI-HC-RFP-20	20	4°C	Cycle	1 Swing Solid	5 Adj.	30 ¾	29	74 5/8
 DAI-HC-SLB-23	23	4°C	Cycle	1 Swing Solid	4 Adj.	26 7/8	34 ¾	81 ¾*
 DAI-HC-SLB-26	26	4°C	Cycle	1 Swing Solid	4 Adj.	28 ¾	36 ½	81 ¾*
 DAI-HC-SLB-36	36	4°C	Cycle	2 Swing Solid	8 Adj.	39 5/8	34 ¾	81 ¾*
 DAI-HC-SLB-49	49	4°C	Cycle	2 Swing Solid	8 Adj.	54	34 ¾	81 ¾*
 DAI-HC-SLB-72	72	4°C	Cycle	3 Swing Solid	12 Adj.	81	34	81 ¾*

Specifications are subject to change

*Height includes 3 ½" casters