

The General Purpose Auto Defrost Freezers are intended for general purpose applications where precise temperature control and maintenance are not essential and there are no requirements for temperature alarms.

†DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assists in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer. These air temperature spikes can be harmful to some medical and laboratory products. Auto defrost units have programmed defrost cycles that will reach air temperatures of approximately 0°C.



DAI-HC-AFB-1420

DAI-HC-AFB-1720

FEATURING

- Mechanical thermostat
- Adjustable temperature range: -15°C to -20°C
- Auto defrost
- Remote alarm contacts
- Keyed door lock
- Magnetic door gasket for positive seal
- Recessed door handle
- Leveling legs
- Probe access port
- One year parts and labor warranty, plus an additional four year compressor parts warranty



Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates.

Freezer Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
†DAI-HC-AFB-1420	14	-20°C	Auto	1 Solid	4 Adj. + Basket	28	26 ¼	67 ¾
†DAI-HC-AFB-1720	17	-20°C	Auto	1 Solid	4 Adj. + Basket	28	29 ¾	68

Specifications are subject to change

†DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assists in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.