

General Purpose Flammable Storage Refrigerators & Freezer meet NFPA and OSHA guidelines 45 and 70 for the refrigerated storage of flammable (volatile) materials in general laboratory work areas. There are no electrical components located in the interior of the unit and the compressor's electrical components have been sealed in a vapor-proof enclosure for additional safety. **These units are NOT designed for use in environments where volatile/explosive conditions could potentially exist.**

Freezers feature a manual defrost operating system to limit product-damaging temperature spikes. **Periodic manual defrosting is required.**

FEATURING

- Mechanical thermostat
- Units meet NFPA guidelines 45 and 70
- Adjustable operating range:
Refrigerator: 1°C to 10°C
Freezer: -15°C to -25°C
- Defrost:
Refrigerator: Cycle defrost
Freezer: Manual defrost
- Keyed door lock
- Magnetic door gasket for positive seal
- Heavy duty door handle
- Leveling legs
- Probe access port
- CFC-free refrigerant, foam insulation and packaging
- UL/C-UL listed
- One year parts and labor warranty



DAI-FFB-30



DAI-FRCB-30

Refrigerator Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
DAI-FRB-30	30	4°C	Cycle	1 Solid	5 Adj.	35 ¾	35 ½	72 ½
Combination Model	Ft ³ (Ref/Frz)	Set Point (Ref/Frz)	Defrost (Ref/Frz)	Doors (Ref/Frz)	Shelves (Ref/Frz)	Width"	Depth"	Height"
DAI-FRCB-30	30	4°C/-20°C	Cycle/Manual	1 Solid/N/A	3 Adj./N/A	35 ¾	35 ½	72 ½
Freezer Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
DAI-FFB-30	30	-20°C	Manual	1 Solid	4 Fixed	35 ¾	35 ½	72 ½

Specifications are subject to change