

General Purpose Freestanding Undercounter/Countertop Refrigerator and Freezers

General Purpose Laboratory Freestanding Refrigerator and Freezers are intended for General Purpose Applications where the precision of a microprocessor temperature controller and security alarms are not required. **For freezers, periodic manual defrosting is required.**

General Purpose products are not intended for high value laboratory supplies, biologicals and pharmaceuticals where precise temperature controls, maintenance and monitoring displays/alarms are required.



DAI-HC-UCFS-0504W




DAI-HC-UCFS-0220M

FEATURING

- Mechanical thermostat
- Adjustable Temperature Range:
Refrigerator: 1°C to 10°C
Freezer: -15°C to -25°C
- Defrost:
Refrigerator: Cycle
Freezer: Manual
- Forced draft air circulation (Refrigerator only)
- Keyed door lock
- Magnetic door gasket for positive seal
- Probe access port
- UL/C-UL listed
- One year parts and labor warranty
- These models require 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure



Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates.

Refrigerator Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 DAI-HC-UCFS-0504W	5	4°C	Cycle	1 Solid	2 Adj./ 1 Fixed	20 ¼	22	31 ½
Freezer Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 DAI-HC-UCFS-0220M	1.5	-20°C	Manual	1 Solid	1 Fixed	18 ½	19 ½	19 ½
 DAI-HC-UCFS-0420W	4	-20°C	Manual	1 Solid	2 Fixed	20 ¼	21 ⅝	31 ½

Specifications are subject to change