

# Pharmacy Plus Series

## Controlled Auto Defrost Freezers and Combinations

Our proprietary "Controlled Auto Defrost" technology provides the convenience of an auto-defrost design while still ensuring that vaccines remain at consistent and stable temperatures below 5°F (-15°C).

Convenient frost removal in freezer while maintaining "spike-free" temperature ranges. No need to relocate samples as would be required during a manual defrost process. No worry of product-damaging temperature excursions as would be experienced with an old technology/traditional auto defrost freezer.



PH-DAI-HC-RFC7A-CAD



PH-DAI-HC-UCFS-0120A-CAD








- Digital temperature display
- Audible and visual temperature alarms
- Resettable Min/Max
- °C/°F switchable



Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates.

### FEATURING

- Microprocessor temperature controller
- Adjustable temperature range:  
Refrigerator: 36°F to 46°F (2°C to 8°C)  
Freezer (-4°F): 4°F to -13°F (-15°C to -25°C)
- Remote alarm contacts, refrigerator and freezer
- Back wall evaporator cover (refrigerators only)
- Forced draft circulation (refrigerators only)
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- LED interior lights
- Probe access ports, refrigerator and freezer
- Independent mechanical refrigeration systems with separate compressors, refrigerator and freezer
- Temperature Monitor Device - complies with the current CDC guidelines, 3 years certification of calibration, "buffered" probe in the product simulated solution, min/max memory, °F/°C switchable, field installable, and visual & audible temperature alarms
- Vaccine storage power cord warning label
- Pharmacy refrigerator/freezer toolkit with temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- UL/C-UL listed
- Two year parts and labor warranty, plus an additional three year compressor parts warranty
- These models require 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure

Freezer Model	Ft³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 PH-DAI-HC-UCFS-0120A-CAD	1.7	-4°F	Controlled Auto	1 Solid	2 Adj.	23 ¾	24	20 ⅜
Combination Model	Ft³ (Ref/Frz)	Set Point (Ref/Frz)	Defrost (Ref/Frz)	Doors (Ref/Frz)	Shelves (Ref/Frz)	Width"	Depth"	Height"
 PH-DAI-HC-RFC7SA-CAD	5.2/1	41°F/-4°F	Cycle/Controlled Auto	2 Solid	3 Adj./2 Adj.	23 ¾	24	53 ⅛
 PH-DAI-HC-RFC7A-CAD	5.2/1	41°F/-4°F	Cycle/Controlled Auto	Glass/Solid	3 Adj./2 Adj.	23 ⅝	24 ¼	53 ⅛
 PH-DAI-HC-RFC12A-CAD	10.5/1	41°F/-4°F	Cycle/Controlled Auto	2 Solid	6 Adj./2 Adj.	23 ¾	26 ½	80 ¼
 PH-DAI-HC-RFC12GA-CAD	10.5/1	41°F/-4°F	Cycle/Controlled Auto	Glass/Solid	6 Adj./2 Adj.	23 ¾	26 ½	80 ¼

Specifications are subject to change  
\*Left hinge model available