

The Pharmacy Plus Series Freestanding Undercounter/Countertop Freezers feature a microprocessor temperature controller and display, which allows for precise temperature control, verification and recovery versus mechanical style thermostats. **For manual defrost units, periodic manual defrosting is required.**




PH-DAI-HC-UCFS-0520

PH-DAI-HC-UCFS-0120








- Digital temperature display
- Audible and visual temperature alarms
- Resettable Min/Max
- °C/°F switchable

 Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates.

### FEATURING

- Microprocessor temperature controller
- Adjustable temperature range:  
(-4°F): 4°F to -13°F (-15°C to -25°C)  
(-22°F): -13°F to -31°F (-25°C to -35°C)  
(-40°F): -31°F to -49°F (-35°C to -45°C)
- Auto & Manual defrost available
- Remote alarm contacts
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- Probe access port
- Temperature Monitor Device - complies with the current CDC guidelines, 3 years certification of calibration, "buffered" probe in the product simulated solution, min/max memory, °F/°C switchable, field installable, and visual & audible temperature alarms
- Vaccine storage power cord warning label
- Pharmacy refrigerator/freezer toolkit with temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- UL/C-UL listed
- Two year parts and labor warranty, plus an additional three year compressor parts warranty
- These models require 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure

Freezer Model	Ft <sup>3</sup>	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 †PH-DAI-HC-UCFS-0120A*	1.3	-4°F	Auto	1 Solid	2 Adj.	23 ¾	24	20 ¾
 PH-DAI-HC-UCFS-0120*	1.7	-4°F	Manual	1 Solid	2 Adj.	23 ¾	24	20 ¾
 PH-DAI-HC-UCFS-0520	4	-4°F	Manual	1 Solid	2 Fixed + 1 Basket	23 ¾	24 ¼	32 ¾
 PH-DAI-HC-UCFS-0430	4	-22°F	Manual	1 Solid	3 Adj.	23 ¾	24	32 ½
 PH-DAI-HC-UCFS-0440	4	-40°F	Manual	1 Solid	3 Adj.	23 ¾	24	32 ½

Specifications are subject to change

\*Left hinge model available

†DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assists in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.