


The Plus Series Built-in Dual Temperature Refrigerator/Freezers Combinations feature separate microprocessor temperature controllers allowing for precise temperature control, verification and recovery. Refrigerators feature an interior fan driven, forced-draft air circulation system provides excellent temperature uniformity and recovery after door openings. **For manual defrost units, periodic manual defrosting is required.**

Models are designed with a “front-breathing” ventilation system, eliminating the need for clearance around the top and sides of the unit; enabling the user to incorporate these units into cabinets and other applications where free space is limited.



DAI-HC-RFC9SS

DAI-HC-RFC9G






Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates.

FEATURING

- Built-in, front breathing technology
- Microprocessor temperature controller, refrigerator and freezer
- Adjustable temperature range:
Refrigerator: 1°C to 10°C
Freezer: -15°C to -25°C
- Defrost:
Refrigerator: Cycle
Freezer: Manual
- Remote alarm contacts, refrigerator and freezer
- Back wall evaporator cover (refrigerators only)
- Forced draft circulation (refrigerators only)
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- Independent mechanical refrigeration systems with separate compressors
- Two probe access ports, refrigerator and freezer
- UL/C-UL listed
- Two year parts and labor warranty, plus an additional three year compressor parts warranty



- Digital temperature display
- Audible and visual temperature alarms
- Resettable Min/Max
- °C/°F switchable

Combination Model	Ft ³ (Ref/Frz)	Set Point (Ref/Frz)	Defrost (Ref/Frz)	Doors (Ref/Frz)	Shelves (Ref/Frz)	Width"	Depth"	Height"
 DAI-HC-RFC9SS*	4.5/4.2	4°C/-20°C	Cycle/Manual	2 Solid Stainless	2 Adj., 1 Fixed/2 Fixed	23 ¾	24 ½	66 ¾
 DAI-HC-RFC9*	4.6/4.2	4°C/-20°C	Cycle/Manual	2 Solid	2 Adj., 1 Fixed/2 Fixed	23 ¾	24 ½	66 ¾
 DAI-HC-RFC9G*	4.6/4.2	4°C/-20°C	Cycle/Manual	Glass/Solid	2 Adj., 1 Fixed/2 Fixed	23 ¾	25 ½	66 ¾

Specifications are subject to change
*Left hinge model available