

The flammable material storage refrigerator and freezer feature a digital display, alarms and remote alarm contacts. These units meet NFPA and OSHA guidelines 45 and 70 for the refrigerated storage of flammable (volatile) materials in general laboratory work areas.

There are no electrical components located in the interior of the unit and the compressor's electrical components have been sealed in a vapor-proof enclosure for additional safety. **These units are NOT designed for use in environments where volatile / explosive conditions could potentially exist.**

The freezer features a manual defrost operating system to limit product-damaging temperature spikes found in auto defrost units. **Periodic manual defrosting is required.**

FEATURING

- Microprocessor temperature controller
- Units meet NFPA guidelines 45 and 70
- Adjustable temperature range:
Refrigerator: 1°C to 10°C
Freezer: -15°C to -25°C
- Defrost:
Refrigerator: Cycle defrost
Freezer: Manual defrost
- Remote alarm contacts
- Keyed door lock
- Magnetic door gasket for positive seal
- Heavy duty door handle
- Leveling legs
- Probe access port
- CFC-free refrigerant, foam insulation and packaging
- UL/C-UL listed
- Two year parts and labor warranty



DAI-FFP-30

DAI-FRP-30



- Digital temperature display and alarm module with battery back-up to support temperature monitoring and alarms
- Sample simulated digital temperature display; .1° Resolution, °C/°F switchable
- Full suite of alarms including - audible and visual high/low temperature alarms, sensor error alarm, power failure alarm
- Resettable 1-point min/max history
- Alarm validation (audible, visual, and remote alarm)
- Alarm mute with ring back

Refrigerator Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
DAI-FRP-30	30	4°C	Cycle	1 Solid	5 Adj.	35 ¾	35	78 ½
Freezer Model								
DAI-FFP-30	30	-20°C	Manual	1 Solid	4 Fixed	35 ¾	35	78 ½

Specifications are subject to change