

The Plus Series Freestanding Undercounter/Countertop Refrigerators feature a microprocessor temperature controller and display, which allows for precise temperature control, verification and recovery versus mechanical style thermostats. Refrigerators feature an interior fan driven, forced-draft air circulation system provides excellent temperature uniformity and recovery after door openings.




DAI-HC-UCFS-0504



DAI-HC-UCFS-0104G

## FEATURING

- Microprocessor temperature controller
- Adjustable temperature range: Refrigerator: 1°C to 10°C
- Cycle defrost
- Remote alarm contacts
- Back wall evaporator cover
- Forced draft circulation
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- LED interior lights (select models)
- Probe access port
- UL/C-UL listed
- Two year parts and labor warranty, plus an additional three year compressor parts warranty
- These models require 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure



- Digital temperature display
- Audible and visual temperature alarms
- Resettable Min/Max
- °C/°F switchable

Glass Door Model	Ft <sup>3</sup>	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 DAI-HC-UCFS-0104*	1	4°C	Cycle	1 Solid	2 Adj.	17 ¼	19 ¼	21 ¼
 DAI-HC-UCFS-0104G*	1	4°C	Cycle	1 Glass	2 Adj.	17 ¼	19 ¼	21 ¼
 DAI-HC-UCFS-0204	2.5	4°C	Cycle	1 Solid	4 Adj.	17 ¾	19 ½	28 ⅞
 DAI-HC-UCFS-0204G	2.5	4°C	Cycle	1 Glass	4 Adj.	17 ¾	19 ½	28 ⅞
 DAI-HC-UCFS-0504	5.2	4°C	Cycle	1 Solid	3 Adj.	23 ¾	24	32 ⅞
 DAI-HC-UCFS-0504G	5.2	4°C	Cycle	1 Glass	3 Adj.	23 ¾	24	32 ⅞

Specifications are subject to change  
\*Left hinge model available