

The TempLog Plus Chromatography Refrigerators are purpose built for scientific and medical applications. Units feature a digital microprocessor temperature controller allowing for the precise temperature management necessary for critical samples and supplies. A forced air directional refrigeration system provides superior temperature uniformity and recovery after door openings. An automatic cycle-defrost design prevents product damaging temperature spikes.

The TempLog Touch Screen temperature display and alarm module offers the latest in advanced temperature control and data logging. Industry-leading product security with a monitoring and alarm system that is 100% independent of the temperature controller.

All TempLog Plus Chromatography Refrigerators require two electrical outlets.

FEATURING

- Microprocessor temperature controller
- Adjustable temperature range: 1°C to 10°C
- Cycle defrost
- One GFCI protected duplex electrical outlet
- Two 2" access ports with spring loaded self-closing, with gasketed metal covers
- Remote alarm contacts
- Forced draft air circulation
- Keyed door lock
- Magnetic door gasket for positive seal
- Heavy duty door handle
- LED interior lights
- Swivel casters (3 ½")
- Probe access port
- ETL Listed
- Three year parts and labor warranty, plus an additional two year compressor parts warranty



DAI-HC-CP-47-TS

DAI-HC-CP-23-TS





DAI-HC-CP-69-TS



- Password-protected touchscreen digital temperature display and alarm module with battery back-up to support temperature monitoring and alarms
- Sample simulated digital temperature display; .1° Resolution, °C/°F switchable
- Active temperature scale bar display
- Resettable 1-point min/max history
- Digital 24-hour temperature chart
- Temperature data logging
- USB port for data transfer and software updates
- Real-time operational display clock
- 4-20mA output
- Full suite of alarms with validation including: audible and visual high/low temperature alarms, sensor error alarm, power failure alarm, door ajar alarm, low battery alarm
- Operation indicator
- Alarm mute with ring back



Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates.

Refrigerator Model	Ft ³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 DAI-HC-CP-16-TS	16	4°C	Cycle	1 Swing Glass	3 Adj.	25 ¾	29 ¾	79*
 DAI-HC-CP-23-TS	23	4°C	Cycle	1 Swing Glass	4 Adj.	27 ⅝	34 ¾	81 ¾*
 DAI-HC-CP-26-TS	26	4°C	Cycle	1 Swing Glass	4 Adj.	29 ⅞	36 ½	81 ¾*
 DAI-HC-CP-33-TS	33	4°C	Cycle	2 Sliding Glass	8 Adj.	40 ⅜	33 ¼	81 ¾*
 DAI-HC-CP-36-TS	36	4°C	Cycle	2 Swing Glass	8 Adj.	40 ⅜	34 ¾	81 ¾*
 DAI-HC-CP-47-TS	47	4°C	Cycle	2 Sliding Glass	8 Adj.	54 ⅞	33 ¼	81 ¾*
 DAI-HC-CP-49-TS	49	4°C	Cycle	2 Swing Glass	8 Adj.	54 ¾	34 ¾	81 ¾*
 DAI-HC-CP-69-TS	69	4°C	Cycle	3 Sliding Glass	12 Adj.	81 ¾	33 ¼	81 ¾*
 DAI-HC-CP-72-TS	72	4°C	Cycle	3 Swing Glass	12 Adj.	81 ¾	34	81 ¾*

Specifications are subject to change

*Height includes 3 ½" casters