


Accessories

1400 LXP



Basic basket with reversible tray lid: Catalog No. PST and GC



Basic basket with spray arm: Catalog No. PSBT






Basket for culture dishes: Catalog No. BP




Baskets	1400 LXP	Part number
Useful for beakers, trays, dishes, dissolution vessels, funnels, instruments and jars.		
Basic basket	PST 14	02600051
Basic basket with spray arm	PSBT 14	02600050
Basket, half with spray arm	PSBT 1/2 14	02600048
Support for culture plates (20 pcs./support max 3 supports) to be used with the PST basket	20 PL	01030417

Basket covers for basic grid baskets	1400 LXP	Part number
Accessory basket covers are available to secure small or lightweight labware during the wash cycle.		
Reversible tray lid for PST/PSBT	GC 14	02250017
Reversible tray lid, half size for PST/PSBT	GC 1/2 14	02250024



Short injectors	1400 LXP	Part number
For inside out cleaning of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.		
48 short jets, 6 mm diameter	48 IXC 14	02600059
64 short jets, 6 mm diameter	64 IXC 14	02600063
64 short jets, mixed 2 mm and 6 mm diameter	64 IXA/2 14	02600062
81 short jets, 6 mm diameter	81 IXC 14	02600068
81 short jets, 6 mm diameter	81 IXC 14	02600068
81 short jets, mixed 2 mm and 6 mm diameter	81 IXA/2 14	02600066
100 short jets, 2 mm diameter	100 IXA 14	02600054
100 short jets, 6 mm diameter	100 IXC 14	02600053
100 short jets, mixed 2 mm and 6 mm diameter	100 IXA/2 14	02600052





Long injectors	1400 LXP	Part number
For inside out cleaning of large narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.		
12 long jets, 8 mm diameter	12 IXL 14	02600044
24 long jets, 8 mm diameter	24 IXL 14	AA0260016
36 long jets, 6 mm diameter	36 IXL 14	02600057
48 long jets, 6 mm diameter	48 IXL 14	02600060
64 long jets, 6 mm diameter	64 IXL 14	02600064





Catalog No. 9 IXL 14





Catalog No. 12 IXL 14





Catalog No. 24 IXL 14



Catalog No. 36 IXL 14




Catalog No. 48 IXL 14



Catalog No. 64 IXL 14

Injector rack covers		
Injector rack covers are available to secure small or lightweight labware during the wash cycle.	1400 LXP	Part number
Injector rack cover full-size	GCI 14	02250023
Injector rack cover half-size	GCI 1/2 14	02250024
Reversible tray lid for injector rack	GCI 14	02250023
Reversible tray lid, half size for injector rack	GCI 1/2 PS	02250025
Cover for 1/2 injector basket	GCI 1/2 14	02250024
Adjustable cross grid	GA 14	02250051




Catalog No. 1 IXL 14

Catalog No. 2 IXL 14

Catalog No. 4 IXL 14

Extra long injectors, 305 mm (12") tall		
For inside out cleaning of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and large volumetric flasks.	1400 LXP	Part number
1 long jet, 12.7 mm diameter, with spray ball	1 IXL 14	02600043
2 long jets, 12.7 mm diameter, with spray ball	2 IXL 14	02600045
4 long jets, 12.7 mm diameter, with spray ball	4 IXL 14	02600046

Combination short/long injectors		
For inside out cleaning of mixed loads of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.	1400 LXP	Part number
36 short and long jets, 6 mm diameter	36 IXLC 14	02600058
48 short and long jets, 6 mm diameter	48 IXLC 14	02600061
64 short and long jets, 6 mm diameter	64 IXLC 14	02600065
81 short and long jets, 6 mm diameter	81 IXLC 14	02600069



Catalog No. 36 IXLC 14

Catalog No. 48 IXL 14

Catalog No. 64 IXLC 14

