


Accessories

910 LX / 1300 LX



Basic basket with reversible tray lid:
Catalog No. PST and GC



Basic basket with spray arm:
Catalog No. PSBT




Basket for culture dishes: Catalog No. BP

Baskets		
Useful for beakers, trays, dishes, dissolution vessels, funnels, instruments and jars.		
	910 LX / 1300 LX	Part number
Basic basket	PST	02600011
Basic basket with spray arm	PSBT	02600010
Support for culture plates (20 pcs./support max 3 supports) to be used with the PST basket	20 PL	01030417

Basket covers		
Accessory basket covers are available to secure small or lightweight labware during the wash cycle.		
	910 LX / 1300 LX	Part number
Reversible tray lid for PST/PSBT	GC	01030569
Reversible tray lid, half size for PST/PSBT	GC1/2	01030570
Reversible tray lid for injector rack	GC	01030585
Reversible tray lid, half size for injector rack	GC1/2	01030586




Short injectors		
For inside out cleaning of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.	910 LX / 1300 LX	Part number
Half rack, 27 short jets, 6 mm diameter	27 IXC 1/2	02600022
36 short jets, 6 mm diameter	36 IXC	02600023
41 short jets, 6 mm BOD rack	41 IXC	02600003
48 short jets, 2 mm diameter	48 IXA	02600026
48 short jets, 6 mm diameter	48 IXC	02600027
64 short jets, 2 mm diameter	64 IXA	02600030
64 short jets, 6 mm diameter	64 IXC	02600032
64 short jets, mixed 2 mm and 6 mm diameter	64 IXA/2	02600031
100 short jets, 2 mm diameter	100 IXA	02600016
100 short jets, 6 mm diameter	100 IXC	02600018
100 short jets, mixed 2 mm and 6 mm diameter	100 IXA/2	02600017




Catalog No. 36 IXC

Catalog No. 48 IXC; IXA


Catalog No. 64 IXA; IXC; IXA/2



Catalog No. 100 IXA; IXC; IXA/2




IXC: Injector rack, short injectors, 3" (76 mm) tall, jets 6 mm diameter




IXA: Injector rack, short injectors, 3" (76 mm) tall, jets 2 mm diameter

Long injectors	910 LX / 1300 LX	Part number
For inside out cleaning of large narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.		
5 long jets, 12.7 mm diameter, with spray ball	5 IXL	02600007
9 long jets, 12.7 mm diameter, with spray ball	9 IXL	AA02250223
12 long jets, 8 mm diameter	12 IXL	AA02250224
20 long jets, 8 mm diameter	20 IXL	02600021
36 long jets, 6 mm diameter	36 IXL	02600024
48 long jets, 6 mm diameter	48 IXL	02600028
64 long jets, 6 mm diameter	64 IXL	02600033
100 long jets, 6 mm diameter	100 IXL	02600019




Catalog No. 5 IXL




Catalog No. 9 IXL




Adjustable clip positions labware on injector for optimal cleaning performance.




Catalog No. 12 IXL




Catalog No. 20 IXL




Catalog No. 36 IXL



Catalog No. 48 IXL




Catalog No. 64 IXL




Catalog No. 100 IXL


Extra long injectors For inside out cleaning of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and large volumetric flasks.		
	910 LX / 1300 LX	Part number
1 long jet, 12.7 mm diameter, with spray ball	1 IXL	02600001
2 long jets, 12.7 mm diameter, with spray ball	2 IXL	02600014
4 long jets, 12.7 mm diameter, with spray ball	4 IXL	02600005




Catalog No. 1 IXL




Catalog No. 2 IXL




Catalog No. 4 IXL









Catalog No. 36 IXLC



Catalog No. 48 IXLC; IXCLA




Catalog No. 64 IXLC



Combination short/long injectors For inside out cleaning of mixed loads of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.		
	910 LX / 1300 LX	Part number
36 short and long jets, 6 mm diameter	36 IXLC	02600025
48 short and long jets, 6 mm diameter	48 IXLC	02600029
48 short and long jets, 6 mm diameter, 2 rows 2 mm diameter	48 IXLCA	02600006
64 short and long jets, 6 mm diameter	64 IXLC	02600034
100 short and long jets, 6 mm diameter	100 IXLC	02600020


Injector KIT packages KIT enable customization of standard injector racks to fit specific needs. Each KIT include 10 complete injectors with their supports and nozzles.		Part number
Long injectors, length of 230 mm and 6 mm diameter		6007500984
Long injectors, length of 171 mm and 6 mm in diameter		6007500890
Short injectors, length of 86 mm and 6 mm in diameter		6007500957
Short injectors, length of 56 mm and 6 mm in diameter		6007501012
Long injectors, length of 171 mm and 8 mm in diameter		6007501011
Needle injectors, length of 86 mm		6007501007




Multipurpose racks	910 LX / 1300 LX	Part number
For a combination of labware and pipettes.		
Basic basket plus 27 short and long jets	PS/27 IXLC	02600039
25 pipettes, angled plus 24 short and long jets	25 PIP/24 IXLC	02600004
50 pipettes, angled plus 34 short and long jets	50 PIP/34 IXLC	02600037
14 short jets plus 4 long injectors angled and perforated	4 EP/14 IXC	02600041
28 short and long jets plus 4 long injectors, angled and perforated	4 EP/28 IXLC	02600036
One third PS plus 21 short and long jets, plus 25 pipettes	PS/21 IXLC/25 PIP	02600038

Bioreactor rack	910 LX / 1300 LX	Part number
Rack for 0.26-0.79 gal (1 -3 L) Bioreactors	BRWR	AA02600113


Pipette racks	910 LX / 1300 LX	Part number
Accessory pipette racks are configured to accommodate reusable pipettes.		
Grommet style, individually loaded, 100 position	100 PIP	02600040



Catalog No. BRWR



Catalog No. PS/27/IXLC



Catalog No. 100 PIP

Storage stand	910 LX / 1300 LX	Part number
An accessory storage stand offers a convenient place to store accessory racks and baskets when not in use.		
Holds up to 5 racks	RP	01040060

Trolley	910 LX	1300 LX
The stainless steel loading trolley transfers the load and mates with the washer to simplify loading and unloading. The trolley includes a retainer lip to keep the floor dry from residual water. Includes locking casters.		
Loading trolley	TN 910 LX	TN 13

