




910 LX 1300 LX 1400 LXP 1600 LXP 1600 LXP/HE

Accessories

910 LX / 1300 LX





Basic basket with reversible tray lid:
Catalog No. PST and GC



Basic basket with spray arm:
Catalog No. PSBT



Basket for culture dishes: Catalog No. BP

Baskets

Useful for beakers, trays, dishes, dissolution vessels, funnels, instruments and jars.

Basic basket

910 LX / 1300 LX

Part number

PST

02600011

Basic basket with spray arm

PSBT

02600010

Support for culture plates (20 pcs./support max 3 supports) to be used with the PST basket

20 PL

01030417

Basket covers

Accessory basket covers are available to secure small or lightweight labware during the wash cycle.

910 LX / 1300 LX

Part number

Reversible tray lid for PST/PSBT

GC

01030569

Reversible tray lid, half size for PST/PSBT

GC1/2

01030570

Reversible tray lid for injector rack

GC

01030585

Reversible tray lid, half size for injector rack



GC1/2

01030586





Short injectors		
For inside out cleaning of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.	910 LX / 1300 LX	Part number
Half rack, 27 short jets, 6 mm diameter	27 IXC 1/2	02600022
36 short jets, 6 mm diameter	36 IXC	02600023
41 short jets, 6 mm BOD rack	41 IXC	02600003
48 short jets, 2 mm diameter	48 IXA	02600026
48 short jets, 6 mm diameter	48 IXC	02600027
64 short jets, 2 mm diameter	64 IXA	02600030
64 short jets, 6 mm diameter	64 IXC	02600032
64 short jets, mixed 2 mm and 6 mm diameter	64 IXA/2	02600031
100 short jets, 2 mm diameter	100 IXA	02600016
100 short jets, 6 mm diameter	100 IXC	02600018
100 short jets, mixed 2 mm and 6 mm diameter	100 IXA/2	02600017









Catalog No. 36 IXC




Catalog No. 48 IXC; IXA




Catalog No. 64 IXA; IXC; IXA/2



Catalog No. 100 IXA; IXC; IXA/2




IXC: Injector rack, short injectors, 3" (76 mm) tall, jets 6 mm diameter



IXA: Injector rack, short injectors, 3" (76 mm) tall, jets 2 mm diameter

Long injectors For inside out cleaning of large narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.		
	910 LX / 1300 LX	Part number
5 long jets, 12.7 mm diameter, with spray ball	5 IXL	02600007
9 long jets, 12.7 mm diameter, with spray ball	9 IXL	AA02250223
12 long jets, 8 mm diameter	12 IXL	AA02250224
20 long jets, 8 mm diameter	20 IXL	02600021
36 long jets, 6 mm diameter	36 IXL	02600024
48 long jets, 6 mm diameter	48 IXL	02600028
64 long jets, 6 mm diameter	64 IXL	02600033
100 long jets, 6 mm diameter	100 IXL	02600019




Catalog No. 5 IXL




Catalog No. 9 IXL




Adjustable clip positions labware on injector for optimal cleaning performance.




Catalog No. 12 IXL




Catalog No. 20 IXL




Catalog No. 36 IXL



Catalog No. 48 IXL



Catalog No. 64 IXL



Catalog No. 100 IXL

Extra long injectors

For inside out cleaning of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and large volumetric flasks.

910 LX / 1300 LX

Part number

1 long jet, 12.7 mm diameter, with spray ball

1 IXL

02600001

2 long jets, 12.7 mm diameter, with spray ball



2 IXL

02600014



4 long jets, 12.7 mm diameter, with spray ball

4 IXL



02600005



Catalog No. 1 IXL



Catalog No. 2 IXL



Catalog No. 4 IXL

Combination short/long injectors

For inside out cleaning of mixed loads of narrow neck glassware, bottles, cylinders, boiling flasks, Erlenmeyer flasks and volumetric flasks.

910 LX / 1300 LX

Part number

36 short and long jets, 6 mm diameter

36 IXLC

02600025

48 short and long jets, 6 mm diameter

48 IXLC

02600029

48 short and long jets, 6 mm diameter, 2 rows 2 mm diameter

48 IXLCA

02600006

64 short and long jets, 6 mm diameter

64 IXLC

02600034

100 short and long jets, 6 mm diameter

100 IXLC

02600020

Injector KIT packages

KIT enable customization of standard injector racks to fit specific needs. Each KIT include 10 complete injectors with their supports and nozzles.

Part number

Long injectors, length of 230 mm and 6 mm diameter

6007500984

Long injectors, length of 171 mm and 6 mm in diameter

6007500890

Short injectors, length of 86 mm and 6 mm in diameter

6007500957

Short injectors, length of 56 mm and 6 mm in diameter





6007501012

Long injectors, length of 171 mm and 8 mm in diameter


6007501011

Needle injectors, length of 86 mm


6007501007



Catalog No. 36 IXLC



Catalog No. 48 IXLC; IXLCA



Catalog No. 64 IXLC

Multipurpose racks For a combination of labware and pipettes.	910 LX / 1300 LX	Part number
Basic basket plus 27 short and long jets	PS/27 IXLC	02600039
25 pipettes, angled plus 24 short and long jets	25 PIP/24 IXLC	02600004
50 pipettes, angled plus 34 short and long jets	50 PIP/34 IXLC	02600037
14 short jets plus 4 long injectors angled and perforated	4 EP/14 IXC	02600041
28 short and long jets plus 4 long injectors, angled and perforated	4 EP/28 IXLC	02600036
One third PS plus 21 short and long jets, plus 25 pipettes	PS/21 IXLC/25 PIP	02600038


Bioreactor rack	910 LX / 1300 LX	Part number
Rack for 0.26-0.79 gal (1-3 L) Bioreactors	BRWR	AA02600113

Pipette racks Accessory pipette racks are configured to accommodate reusable pipettes.	910 LX / 1300 LX	Part number
--	-------------------------	--------------------


Grommet style, individually loaded, 100 position

100 PIP



02600040



Catalog No. BRWR



Catalog No. PS/27/IXLC



Catalog No. 100 PIP

Storage stand

An accessory storage stand offers a convenient place to store accessory racks and baskets when not in use.

910 LX / 1300 LX

Part number

Holds up to 5 racks

RP

01040060

Trolley

The stainless steel loading trolley transfers the load and mates with the washer to simplify loading and unloading. The trolley includes a retainer lip to keep the floor dry from residual water. Includes locking casters.

910 LX

1300 LX

Loading trolley

TN 910 LX

TN 13

