

General Purpose Laboratory products are intended for general purpose applications where the precision of a microprocessor temperature controller and security alarms are not required. **For freezers, periodic manual defrosting is required.**

General Purpose products are not intended for high value laboratory supplies, biologicals and pharmaceuticals where precise temperature controls, maintenance and monitoring displays/alarms are required.

## FEATURING

- Mechanical thermostat
- Adjustable Operating Temperature Range:  
Refrigerator: 1°C to 10°C  
Freezer: -15°C to -25°C
- Defrost:  
Refrigerator: Cycle and Auto Available  
Freezer: Manual
- Forced draft air circulation (Refrigerators only)
- Keyed door lock
- Magnetic door gasket for positive seal
- Stainless steel door handle
- Probe access port
- High density urethane foam cabinet and door insulation
- CFC-free refrigerant, foam insulation and packaging
- UL/C-UL listed
- Two year parts and labor warranty, plus an additional three year compressor parts warranty



DAI-30R

DAI-RFC-30M

Refrigerator Model	Ft <sup>3</sup>	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
†DAI-30R	30	4°C	Auto	1 Solid	5 Adj.	35 ¾	35	72 ½
Combination Model	Ft <sup>3</sup> (Ref/Frz)	Set Point (Ref/Frz)	Defrost (Ref/Frz)	Doors (Ref/Frz)	Shelves (Ref/Frz)	Width"	Depth"	Height"
DAI-RFC-30M	25/5	4°C/-20°C	Cycle/Manual	1 Solid/N/A	3 Adj. / N/A	35 ¾	35 ½	72 ½

Specifications are subject to change