

Product Name:	ARU-3004
Description:	<p>General Purpose Refrigerator  30 Cu. Ft. Capacity  One solid swing door, right hinged, not reversible  Mechanical thermostat  Temperature range: 0°C to 10°C  Five adjustable shelves  Keyed door lock  Magnetic door gasket for positive seal  Stainless steel door handle  Leveling legs  Probe access port (3/8")  Powder coat paint exterior  High density urethane foam cabinet and door insulation  CFC-free refrigerant (R134a), foam insulation and packaging  Accu-test certified  Two year parts and labor warranty, plus an additional three year  compressor parts warranty  Exterior dimensions: 35 3/4"W x 35" D x 72 1/2"H  Auto defrost  115V, 60Hz, 6.8 amps, 1/3 HP  UL/C-UL listed  Shipping weight: 396 lbs.</p>