



Mopec Premium Biopsy Bags

Ensure unrestricted fluid movement around tissue with Mopec Premium Biopsy Bags. Easy to open and available in two convenient sizes, you'll find our nylon mesh biopsy bags are the optimal solution for protecting specimens during tissue processing.


- ✔ Specimen Control
- no trampoline effect
- ✔ Efficient Packaging
- save on waste and damage reduction
- ✔ Multiple Filament Design
- easy to grab and easier to open




UNDER THE MICROSCOPE: WHY CUSTOMERS PREFER MOPEC BIOPSY BAGS



< Top: Seam Welding >



< Bottom: Weaved Filament >



Mopec Premium Biopsy Bag

Multifilament design ensures consistent seam strength.

Competitors Biopsy Bag

Dualfilament and Monofilament design creates an inconsistent seam strength.

MANUFACTURED WITH SPECIMEN SAFETY IN MIND

Unlike other biopsy bags on the market, Mopec Premium Biopsy Bags are built with a twisted, multiple filament design that creates a textured exterior. This specialty pattern gives the surface of the bags a tacky feel, making them easier to grab. The consistent seam strength on the edges of the bag not only make them easier to open, but helps prevents “the trampoline effect,” therefore reducing specimen loss.



EFFICIENT PACKAGING

While other biopsy bags are packaged loosely, Mopec Premium Biopsy Bags are organized specifically to make it easier for the end user. Mopec Biopsy Bags are stacked and oriented in the same direction to reduce opening time and create a more efficient workspace for any pathologist assistant.