

AniGARD®



The AniGARD® AG303 delivers the same trusted protection and performance as our larger units in a more compact, space-conscious design. At just 3 feet wide, this model is ideal for facilities where space is limited — without compromising ISO Class 4 air cleanliness, operator comfort, or allergen containment. Whether you're working in smaller animal rooms or looking to supplement your existing equipment, the AG303 provides an ergonomic, maneuverable solution for high-demand animal research environments.



Engineered for modern animal research labs, the AniGARD® AG303 delivers top-tier protection and comfort in a streamlined 3-foot footprint.



Learn more at
Bakerco.com

FEATURES

- **Compact 3-Foot Width** – Perfect for space-constrained labs or mobile setups, while still accommodating standard cage sizes and procedures.
- **Momentum Air Curtain** – Creates a high-velocity airflow barrier along the work surface perimeter, preventing room air contaminants from entering while ensuring exceptional containment and easy access for users.
- **Adjustable Height** – Features an electric hydraulic lift system that allows the work surface to be set between 34" and 44", ensuring ergonomic comfort for individual operators and accommodating various tasks.
- **Hinged Viewscreens** – Designed to open wide up to 21.75", these screens make it easy to handle large items, facilitate cage changes, and provide excellent visibility of the workspace during use.
- **Energy Efficiency** – Powered by StediFLOW™ technology, the system uses less energy, generates less heat, and extends the life of HEPA filters, reducing maintenance and operational costs over time.

APPLICATION

The AniGARD® AG303 is designed for cage changing, animal handling, and other routine lab procedures where particulate and allergen control are critical. Its smaller footprint makes it especially suited to high-density research environments or satellite facilities, offering full protection performance with minimal space requirements.

TECHNICAL SPECIFICATION ▼

MODEL NUMBER	AG303
Size	3'
Exterior Dimensions (w x d x h [min to max])	46(1/10)" x 32(1/2)" x [77(1/2)" - 81(1/4)]
Interior Dimensions (w x d x h)	36" x 30" x 22(5/16)"
Work Surface Height [min to max]	34" - 44"
Usable Workspace (w x d)	31(5/16)" x 21(13/16)"
Weight	478 lbs
Shipping Weight	600 lbs
Access Opening (working position)	14"
Access Opening (loading/cleaning position)	21(3/4)"
Operating Amperage	3.7 Amps
Power Consumption	426 Watts
Heat Generation	1,454 BTU/hr
Electrical Requirements*	115 V AC, 15 Amps, 60 Hz
Noise - NSF Method	61 dBA
Lighting - Minimum average at work surface	100 foot candles

RESOURCES

Scan the codes to view our helpful resources to maximize the use of the AniGARD® Clean Bench.



Visit the
AniGARD® VF
Product page



Download
the AniGARD®
Brochure



View the
AniGARD®
Airflow demo