

The Controlled Room Temperature Cabinets are purpose built for medical applications. These units feature a microprocessor temperature controller and display allowing for precise temperature management, verification and recovery. An interior fan driven, forced-draft air circulation system provides excellent temperature uniformity and recovery after door openings. Medications are security backed with an extensive package of audible and visual alarms.

Room temperature storage of medication is defined by the United States Pharmacopoeia as 68°F to 77°F (20°C to 25°C). A controlled temperature environment will ensure approved storage; preventing time wasted on unnecessary expense and product recalls. **All Controlled Room Temperature models are designed with cooling functionality only, they are NOT designed with a heating function to protect against temperatures below 68°F (20°C).**



CRT-DAI-HC-10PG

CRT-DAI-HC-10PS





- Digital temperature display
- Audible and visual temperature alarms
- Resettable Min/Max
- °C/°F switchable



Environmentally friendly, natural hydrocarbon (HC) refrigerants that vastly reduce global-warming potential while meeting new EPA/SNAP mandates.

## FEATURING

- Microprocessor temperature controller
- Adjustable temperature range: 68°F to 77°F (20°C to 25°C)
- Cycle defrost
- Remote alarm contacts
- Forced draft air circulation
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- Recessed door handle
- LED interior lights
- Probe access port
- Temperature Monitor Device - complies with the current CDC guidelines, 3 years certification of calibration, "buffered" probe in the product simulated solution, min/max memory, °F/°C switchable, field installable, and visual & audible temperature alarms
- Vaccine storage power cord warning label
- Pharmacy refrigerator/freezer toolkit with temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- UL/C-UL listed
- Two year parts and labor warranty, plus an additional three year compressor parts warranty
- These models require 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure

Glass Door Model	Ft³	Set Point	Defrost	Doors	Shelves	Width"	Depth"	Height"
 CRT-DAI-HC-10PG	10.5	72°F	Cycle	1 Glass	6 Adj.	23 ¾	26 ½	59
Solid Door Model								
 CRT-DAI-HC-10PS	10.5	72°F	Cycle	1 Solid	6 Adj.	23 ¾	26 ½	59

Specifications are subject to change