

Product Name:	PH-DAI-HC-10PS
Description:	Plus Series Pharmacy/Vaccine Compact Refrigerator 10.5 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Calibrated using a NIST traceable device, three year certificate included Remote alarm contacts Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall Bottom footed shelf to allow air circulation Back wall evaporator cover Vaccine storage power cord warning label Forced draft circulation Keyed door lock Magnetic door gasket for positive seal Leveling legs LED interior light Probe access port (3/8") Powder coat exterior High density urethane foam cabinet insulation HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾"W x 26 ½"D x 59"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> UL/C-UL listed ENERGY STAR CERTIFIED Cycle defrost 115V, 60 Hz, 1/3 HP, 2.3 Amps Energy Consumption 2.42 kWh/day Shipping weight: 183 lbs.