

Product Name:	PH-DAI-HC-SSP-49
Description:	<p>Plus Series Pharmacy/Vaccine Stainless Steel Refrigerator</p> <p>49 Cu. Ft. Capacity</p> <p>Two stainless steel swing doors</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display</p> <p>Audible and visual high and low temperature alarms</p> <p>Remote alarm contacts</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Adjustable operating temperature range: 36°F – 46°F (2°C- 8°C)</p> <p>Calibrated using a NIST traceable device, three year certificate included</p> <p>Six adjustable shelves</p> <p>Vaccine storage power cord warning label</p> <p>Forced draft air circulation</p> <p>Keyed door locks</p> <p>Magnetic door gasket for positive seal</p> <p>Integrated door handles</p> <p>Swivel casters (adds 6 ½" to height)</p> <p>Interior light</p> <p>Probe access port (3/8")</p> <p>Stainless steel interior floor</p> <p>Stainless steel exterior</p> <p>High density urethane foam cabinet insulation</p> <p>Hydrocarbon, natural refrigerant (R290)</p> <p>Pharmacy refrigerator/freezer toolkit and temperature logs</p> <p>Two year parts and labor warranty, plus an additional three year compressor parts warranty</p> <p>Exterior dimensions: 54"W x 35"D x 84"H with casters; 77 ½"H without casters</p> <p>Cycle defrost</p> <p>115V, 60 Hz, 4.5 Amps, 1/4 HP</p> <p>ETL listed</p> <p>Shipping weight: 450 lbs.</p>