

Product Name:	PH-DAI-HC-SSP-49FA
Description:	<p>Plus Series Pharmacy Stainless Steel Freezer 49 Cu. Ft. capacity</p> <p>Two swing solid doors, self-closing, not reversible</p> <p>6 ea. Adjustable Shelves</p> <p>Digital display microprocessor temperature controller</p> <p>Adjustable operating temperature: 5°F – -13°F (-15°C - -25°C)</p> <p>Calibrated using a NIST traceable device, three year certificate included</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Probe attached to evaporator for defrost termination sensing Sensor error alarm</p> <p>Continuous product (simulated) temperature display</p> <p>Remote alarm contacts</p> <p>Forced draft air circulation</p> <p>Keyed door lock</p> <p>Magnetic door gaskets for positive seal</p> <p>Integrated door handle</p> <p>LED interior lights are safety shielded and switch controlled Probe access port (3/4")</p> <p>Stainless steel exterior</p> <p>White-powder coated aluminum interior with 304 stainless steel floor</p> <p>Easily cleaned smooth interior</p> <p>Hydrocarbon natural refrigerant (R290)</p> <p>2 year parts and labor warranty, plus an additional 3 year compressor parts warranty</p> <p>Auto defrost</p> <p>4 swivel casters (5 3/4")</p> <p>115V/60 Hz</p> <p>ETL listed</p> <p>Exterior dimensions: 54" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters</p> <p>115V, 60 Hz, 9 Amps, 3/4 HP</p> <p>ETL listed</p> <p>Shipping weight: 495 lbs.</p> <p>DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.</p>