

Product Name:	PH-DAI-HC-UCBI-0420SS-LH
Description:	Plus Series Pharmacy/Vaccine Undercounter Freezer Built-In-LEFT HINGED 4.2 Cu. Ft. Capacity One stainless steel swing door, left hinged Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 5°F – -13°F (-15°C – -25°C) Calibrated using a NIST traceable device, three year certificate included Two fixed shelves Vaccine storage power cord warning label Keyed door lock Recessed door handle Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Stainless steel exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾"W x 24 ½"D x 33 ⅜"H Manual defrost 115V, 60 Hz, 1.5 Amps, 1/8 HP Energy Consumption 0.97 kWh/day UL/C-UL listed ENERGY STAR CERTIFIED Shipping weight: 131 lbs.