

Product Name:	PH-DAI-HC-UCFS-0430
Description:	<p>Plus Series Undercounter Pharmacy Freezer Freestanding (-30°C)  4 Cu. Ft. Capacity  One solid swing door, right hinged, not reversible  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range:-4°F to -22°F (-20°C to -30°C)  Factory Preset at -30°C  Calibrated using a NIST traceable device, three year certificate included  Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall  Bottom footed shelf to allow air circulation  Back wall evaporator cover  Keyed door lock  Magnetic door gasket for positive seal  Recessed door handle  Leveling legs  LED interior light  Probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R290), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Pyxis®, compressor parts only Omnicell® and warranty AcuDose RX® compatible  Exterior dimensions: 23 3/4" W x 24"D x 32 1/8"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost  115V, 60 Hz, 2.2 Amps, 1/4 HP  UL/C-UL listed  ENERGY STAR CERTIFIED  Shipping weight: 128 lbs.</p>