

Product Name:	PH-DAI-RFC-16A
Description:	<p>Plus Series Pharmacy/Vaccine Refrigerator & Freezer Combination 16 Cu. Ft. capacity (13 cu. ft. refrigerator/3 cu. ft. freezer) Two exterior solid swing, double insulated doors, right hinged, not reversible Individual microprocessor temperature controllers for both the refrigerator and freezer Independent temperature displays Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – 13°C (-15°C - -25°C Calibrated using a NIST traceable device, three year certificate included Audible and visual high and low temperature alarms, refrigerator and freezer Sample temperature in freezer does not rise above -20°C during defrost cycle Uniformity of samples throughout the freezer ($\pm 2^{\circ}\text{C}$ throughout chamber) Uniformity of samples throughout the refrigerator ($\pm 1^{\circ}\text{C}$ throughout chamber) 36 hour battery back-up on alarms, high/low temperature & loss of power Two isolated audible & visual high/low alarm systems remote alarm contacts Remote alarm contacts, refrigerator & freezer 3 adjustable, vented stainless steel shelves refrigerator/ 1 fixed shelf freezer Vaccine storage power cord warning label Leveling legs Keyed door locks, refrigerator & freezer Magnetic door gasket for positive seal Heavy duty door handles Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer) White powder coated exterior CFC-free refrigerant (R134a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 29 ½"W x 34 ¾"D x 78"H Auto defrost freezer / cycle defrost refrigerator 115 V, 60Hz, 7.5 Amps, 1/8 HP UL/C-UL listed Shipping weight: 323 lbs. DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer</p>